en **Instructions for use/Technical description**

USA Caiman 5 Sealing & cutting Instruments

Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusaifus.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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Caiman 5 instruments Sealing & Cutting

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Caiman 5 Sealing & Cutting-instrumenten

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Caiman 5 Sealing & Cutting instrumenter

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Caiman 5 Sealing & Cutting-instrument

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Caiman 5 Sealing & Cutting instrumendid

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Instrumenti za zapiranje in rezanje Caiman 5

hr **Upute za uporabu/Tehnički opis**

Caiman 5 Sealing & Cutting instrumenti

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
Legend

Instrument

- 1 HF activation button
- 2 Rotation wheel
- 3 Moveable jaw part with marking
- 4 Blade activation lever
- 5 Jaw part activation lever
- 6 HF cable with connector

Generator

- 7 See Instructions for Use Lektrafuse HF generator
- 8 Connection socket (foot switch)
- 9 Signal lamp Regrasp
- 10 Display
- 11 Signal lamp Error
- 12 Connection socket (instrument)
- 13 Indicator lamp power-on
- 14 Foot switch

1. About this document

Note

General risk factors associated with surgical procedures are not described in these instructions for use.

1.1 Scope

Aesculap Caiman 5 Sealing & cutting instruments are sterile disposable instruments. In conjunction with the LEKTRAFUSE GENERATOR the Caiman 5 instruments form an HF system for the application of HF energy. Caiman 5 Instruments can only be used in conjunction with the LEKTRAFUSE GENERATOR.

- For article specific instructions for use and material compatibility and lifetime information, see B. Braun eIFU at eifu.bbraun.com

1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

⚠ DANGER

Indicates a possible threat of danger. If not avoided, death or serious injury may result.

⚠ WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

⚠ CAUTION

Indicates a possible threat of material damage. If not avoided, the product may be damaged.

2. Clinical use

2.1 Areas of use and limitations of use

2.1.1 Intended use

Caiman Seal & Cut is a bipolar HF sealing system, which consists of the Lektrafuse HF generator and Caiman instruments.

The Lektrafuse HF generator is used to convert mains power to HF power based on advanced bipolar technology. Caiman Seal & Cut instruments are only compatible with the Lektrafuse HF generator.

The Caiman Seal & Cut system can be used for grasping, preparation, sealing and cutting of tissue, on vessels and vessel bundles with diameters up to and including 7 mm as well as soft tissue.

The Caiman Seal & Cut system can be used during open and minimally invasive surgery in general surgical procedures including, but not limited to surgical specialties such as gynecology, urology, bariatric, colorectal and thoracic surgery.

Refer to the Lektrafuse HF generator's IFU for additional indications, warnings and contraindications.

2.1.2 Indications

Note

The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

For indications, see Intended use.

2.1.3 Absolute contraindications

- Caiman Seal & Cut is not suitable for grasping, preparation, sealing and cutting of hard tissue such as bone or cartilage.
- Caiman Seal & Cut is not approved for tubal sterilization, coagulation of the fallopian tubes, etc.
- Caiman Seal & Cut is not approved for application on parenchyma tissue of organs.

2.1.4 Relative contraindications

Particular caution must be exercised in case of pathological changes of vessels such as atherosclerosis.

Since the safe and effective use of bipolar HF sealing systems strongly depends on factors such as tissue type, pathology and surgical procedure, the above information can only be used as a general guidance. The clinically successful application is dependent upon the knowledge and experience of the surgeon, who is responsible for deciding which structures can be reasonably treated and whether the patient's condition, e.g. the coagulation status, allows a successful treatment while taking into account the safety and warning instructions given in these instructions for use.

2.2 Safety information

2.2.1 Clinical user

General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- Use the product only according to these instructions for use.
- Follow the safety and maintenance instructions.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to use, check that the product is in good working order.
- Keep the instructions for use accessible for the user.

Note

The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

2.2.2 Product specific safety information

- Only combine Aesculap products with each other.
- Always adhere to applicable standards.

Sparks may be created during proper use of the HF instrument, which may lead to ignition or explosion of flammable gases.

- Observe the safety instructions in the instructions for use of the Lektrafuse HF device.

To avoid inadvertent injury and burns to patients and medical personnel, a basic understanding of the surgical HF principles and techniques is required.

To avoid HF burns:

- Always keep the working end of the product in the user's field of vision whenever the HF power is activated.
- Prior to activating the Lektrafuse HF generator, check that the working end of the product is not in contact with any electrically conductive accessories.
- Before each use, inspect products visually for: damage and surface changes on the insulation.
- Never place the product on or next to the patient.
- Observe the safety instructions in the instructions for use of the Lektrafuse HF device.

The rated accessory voltage of the product is 240 Vp.

2.2.3 Sterile packaged products

Sterile products

The product is sterilized using irradiation and sterile packaged.

The safe sterile provision of the product is only guaranteed if the sterile packaging is undamaged and unopened and the use-by date has not passed.

- Store products in their original packaging in a dust-protected, dry and temperature-controlled room.
- Protect products from direct sunlight.
- Remove products from their original protective packaging only just prior to application.
- Visually inspect the product packaging to ensure that the sterile barrier system is intact.
- Do not use products from open or damaged sterile packaging.
- Do not use the product after its use-by date.

2.3 Safe handling and preparation

⚠ DANGER

Risk to the patient and/or user when using the instruments with HF sources other than the Lektrafuse HF generator!

Non-observance can result in injury or death and restriction in the function of the instruments.

- Use the instruments only with the Lektrafuse HF generator.

⚠ DANGER

Failure to observe the instructions below can result in death or injury to the patient!

- Use instruments with caution in patients with pacemakers or defibrillators.

⚠ DANGER

Improper handling may result in injury or damage to the instrument!

- Do not activate the HF current until the instrument has been fully closed and locked.
- Do not open instruments during HF operation, as this will automatically interrupt the sealing.
- Do not activate the blade during the HF process. Insufficient sealing or damage to the instrument may occur.

⚠ WARNING

The use of a damaged or defective instrument may result in injury to the patient and/or user!

- Before using the instrument for the first time, close and lock the jaws and then actuate the blade activation lever to check it functions correctly.
- Check the product for loose, bent, broken, cracked, worn or missing parts before every use.
- Do not use the product if it is damaged or defective.
- If unusual performance or damage to the instrument is observed, discontinue using the instrument immediately. Replace the instrument with a different one of the same type.

⚠ WARNING

The evaporation of tissue fluid can cause injuries!

Tissue fluid can evaporate during sealing because of the HF energy. This may result in collateral damage to the immediately adjacent tissue.

- Aspirate liquids before activating the instrument.

⚠ WARNING

Risk of injury from ignition or explosion of flammable gases! Sparks may occur in the normal use of the HF device.

- Do not use the device in explosion-hazard zones.
- Avoid or aspirate flammable anesthetics and combustible gases (e.g. nitrous oxide, oxygen) during surgery of the head and thorax.
- If possible, use incombustible cleaning and disinfecting agents.
- If flammable cleaning and disinfecting agents or solvents are used (e.g. alcohol, ether): Verify that such agents have evaporated prior to commencing HF surgery.
- Ensure that no combustible fluids accumulate under the patient or in body cavities (e.g. the vagina). Wipe off all fluids before using the HF device.
- Ensure the absence of any endogenous, combustible gases.
- Ensure that any materials soaked in liquid oxygen (e.g. cotton wool, gauze) are kept sufficiently far away from the HF environment to prevent them from igniting.

⚠ WARNING

Dirty electrode surfaces may cause insufficient sealing!

- Keep the instrument's contact surfaces clean during surgery, see chapter Intraoperative cleaning.

⚠ WARNING

The use of incompatible or damaged trocars may result in injury to the patient and/or damage of the instrument!

- When using the instrument through a trocar:
 - Check the compatibility of the trocar and the instrument.
 - Check the trocar for sharp edges and damage.
- Do not use the instrument with damaged trocars.

⚠ WARNING

Removing the instrument through a trocar may result in injury to the patient and/or damage of the instrument!

- Close and lock the instrument before carefully inserting it into the trocar or removing it therefrom.

2.4 Application

DANGER

Risk of injury due to application of the product outside the field of view!

- Use and activate the working end only under visual control.
- Make sure there is no contact with electrically conductive objects before activating the product.
- Do not place the product on or next to the patient.

WARNING

Risk of tissue damage or damage to the instrument due to improper use of the instrument with closed and locked jaws!

- Do not perform the process listed below when the jaws are locked:
 - Shaft rotation
- Make sure that the jaws are open before turning the shaft.

CAUTION

Risk of inadequate sealing or cutting of the tissue!

- Do not bundle or fold the tissue to be sealed, or cut and clamp more than one layer between the jaws.
- Grip the tissue to be sealed centrally between the jaws.
- Do not start the procedure if the jaws of the instrument are in contact with conductive liquid (e.g. blood or saline solution).
- Do not start the procedure if there are conductive objects (e.g. clamps, vessel clamps, clips, etc.) between the jaws of the instrument.

The software of the Lektrafuse HF GENERATOR recognizes the instruments and applies appropriate device settings. The duration of the sealing process is controlled by the generator software.

The sealing process is started by pressing the HF button once and ends automatically after sealing. The sealing process can be interrupted by pressing the HF button again.

Note

In the event of an unsuccessful or incomplete seal and insufficient functionality of the instruments, the generator issues a visual and acoustic alarm.

Note

In the event of visibly inadequate sealing of the tissue, the sealing quality can be improved by reactivating the HF process.

2.4.1 Tissue manipulation, sealing of vessels and tissue bundles and tissue separation

- Remove the instrument from the sterile packaging.
- Insert the connector **6** of the instrument into the connector socket **12** of Lektrafuse the HF generator.
- Use the rotation wheel **2** to move the jaws to the desired position for the intervention. The raised indication on rotation wheel **2** indicates the center position.
- Grip the tissue to be sealed or cut between the jaws.
- Use the jaw part activation lever **5** to close and lock the moveable jaw part **3**.
- To seal the clamped tissue, press the HF activation button **1** at the instrument or the foot switch **14**.

Note

For further options, refer to the instructions for use of the Lektrafuse HF generator.

- Use the blade activation lever **4** to cut the clamped and sealed tissue.

DANGER

Clamping of metallic objects and/or solid, rigid or very thick tissue may cause injury to the patient and/or damage to the instrument!

- Do not grip, seal or cut metal objects such as clips or clamps.
- Do not grip, seal or cut strong, rigid or very thick tissue such as bone or cartilage.
- If an error is reported:
 - Do not cut clamped tissue.
 - Unlock and open the instrument jaws.
 - If the sealing quality is sufficient, the sealed tissue can be cut using surgical scissors or the instrument.
 - If the sealing quality is insufficient, the sealing process must be restarted. If the additional HF cycle results in a normal seal, the tissue can be cut, and the instrument can be opened and removed.
- If repeated error messages occur, discontinue using the instrument immediately. Replace the instrument with a different one of the same type.
- Do not leave the jaw activation lever in the locked position when the instrument is not in use.

2.4.2 Intraoperative cleaning

WARNING

Improper cleaning can cause injury to the patient and/or user!

- Do not activate HF current when cleaning the instrument.

CAUTION

Improper cleaning may result in damage and/or failure of the instrument!

This can result in continuous error messages REGRASP INDICATOR – SHORT (see Instructions for Use Lektrafuse HF generator).

- When cleaning the instrument:
 - Do not immerse the instrument jaws in liquid.
 - Do not use abrasives.
 - Do not use excessive force to clean jaws.
- Clean the instrument only as necessary:
 - If the blade becomes difficult to activate.
 - If the instrument sticks to the tissue.
 - If electrode surfaces become encrusted with tissue residues or body fluids.

When cleaning is required:

- Close and lock the jaws of the instrument and trigger the blade activation lever **4** five to six times. This procedure may loosen tissue residues from the proximal part of the jaws.
- If more thorough cleaning is required, clean jaws with a damp swab and sterile (normal or deionized) water. Remove excess fluid from swab.

Note

Sterile (normal or deionized) water is preferable to saline or chemical cleaning solutions.

After cleaning the instrument:

- Check the inner and outer parts of the jaws for bent or broken parts and other damage.
- Do not use the instrument if it is damaged or defective.
- If repeated error messages occur, immediately discontinue using the instrument. Replace the instrument with a different one of the same type.

3. Single-use products

DANGER

Reuse can cause infections in patients and/or users, and impairment of product functionality. Soiling and/or impaired functions of the product can cause injury, illness or death!

- Do not reprocess the product!

4. Storage

- Store sterile single-use products in germ-proof packaging in a dust-protected, dry, dark and temperature-controlled room.

5. Troubleshooting

Errors, faults and warnings are displayed at the Lektrafuse HF generator as follows:

- An individual tone or a clear tone sequence and interruption of the tone sequence indicating HF operation
- The yellow Regrasp Error light or the red error light on the front of the Lektrafuse HF generator comes on
- Message on the display at the front of Lektrafuse the HF generator
- For detailed information on faults, malfunctions and troubleshooting, refer to the instructions for use of the Lektrafuse HF generator.

6. Maintenance

DANGER

Danger to life of patients and users if the product malfunctions and/or protective measures fail or are not used!

- Do not perform any service or maintenance work while the product is being used on a patient.
- Do not modify the product.

CAUTION

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- Do not modify product.

7. Disposal

WARNING

Risk of infection due to contaminated products!

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

WARNING

Risk of injury due to sharp-edged and/or pointed products!

- When disposing of or recycling the product, ensure that the packaging prevents injury by the product.

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